

# Guide to the Python-based Coil System

<https://github.com/atelier-ritz/CoilSystemPython>

# Table of contents

- Dependencies
- Usage
- Program Structure
- Modify GUI
- Image filters
- Object detection
- Multithreading

# Dependencies


===== Tested on Ubuntu 17.10 =====

- Python 3.6 pre-installed in Ubuntu 17.10
- PyQt5 pip3 install pyqt5
  - What is PyQt <https://riverbankcomputing.com/software/pyqt/intro>
- Opencv pip3 install opencv-python, pip3 install opencv-contrib-python
- Pydc1394
  - Firewire camera module <https://github.com/jordens/pydc1394>
- Qt-designer sudo apt-get install qt4-designer
  - GUI designer


# Usage

Go to the working directory.

Run “python3 main.py”



# Program structure



# Program structure

main.py

callbacks.py *Add your code here*

|  
|  
|

└─syntax.py [highlight the keywords in GUI editor\_vision]

|  
|

└─fieldManager.py [send commands to s826; store XYZ field strength]

|     |     s826.py [control s826 I/O]

|  
|  
|

└─visoin.py [capture frames; apply filters; detect objects]

|     |     filterlib.py [define filters] *Add your code here*


|     |     objectDetection.py [define object detection algorithms] *Add your code here*

|  
|  
|

└─subthread.py [run multithreading tasks] *Add your code here*

# Modify GUI



1. Open "Mainwindow.ui" with qt-designer.



# Modify GUI

1. Open "Mainwindow.ui" with qt-designer.
2. Open "callbacks.py" and edit connectSignals()

More about <signal> and <slot> [http://pyqt.sourceforge.net/Docs/PyQt4/new\\_style\\_signals\\_slots.html](http://pyqt.sourceforge.net/Docs/PyQt4/new_style_signals_slots.html)



<objectName>

```
def connectSignals(self):  
    # General Control Tab  
    self.dsb_x.valueChanged.connect(self.setFieldXYZ)  
    self.dsb_y.valueChanged.connect(self.setFieldXYZ)  
    self.dsb_z.valueChanged.connect(self.setFieldXYZ)  
    self.btn_clearCurrent.clicked.connect(self.clearField)  
    self.dsb_xGradient.valueChanged.connect(self.setFieldXYZGradient)  
    self.dsb_yGradient.valueChanged.connect(self.setFieldXYZGradient)  
    self.dsb_zGradient.valueChanged.connect(self.setFieldXYZGradient)  
    # Vision Tab  
    self.highlighter = syntax.Highlighter(self.editor_vision.document())  
    self.chb_bypassFilters.toggled.connect(self.on_chb_bypassFilters)  
    self.chb_startPauseCapture.toggled.connect(self.on_chb_startPauseCapture)  
    self.btn_refreshFilterRouting.clicked.connect(self.on_btn_refreshFilterRouting)  
    # object detection
```

**self.<objectName>.<signal>.connect(<slot>)**

```
self.chb_startStopSubthread.toggled.connect(self.on_chb_startStopSubthread)  
self.dsb_freq.valueChanged.connect(self.thrd.setFreq)
```



# Image filters

Note:  
ONLY include alphabets, numbers, and  
underbars in your filter name.

1. Open “filterlib.py” and add your custom filter.


Attention: need to handle variable conversion by yourself

E.g. str -> int/float, define upper/lower bounds



2. Use it directly in the GUI.

Double slash to comment it out

Filters are  
connected  
in series  
and applied  
in order.



```
#=====
# canny(minVal,maxVal)
# Input must be a greyscale image
#=====
def canny(inputImage,args):
    arg = args.split(',')
    return cv2.Canny(inputImage,int(arg[0]),int(arg[1]))
```




# Object detection

Note:

When object detection is enabled, the original image overlaid with the detected object will be shown instead of the filtered image.

1. Add the name of your object detection algorithm to the GUI.
2. Define your algorithm in “objectDetection.py”.

*See sample algorithmA() or google “opencv object detection python”*



```
def algorithmA(imageFiltered, imageOriginal, agent):
    nOfSamples = 2
    im2, contours, hierarchy = cv2.findContours(imageFiltered, cv2.RETR_TREE, cv2.CHAIN_APPROX_SIMPLE)
    cnts = sorted(contours, key = cv2.contourArea, reverse = True)[:nOfSamples]
    if len(cnts) > 1:
        targetCnt = cnts[1] # cnt[0] is the edge of the screen
        rect = cv2.minAreaRect(targetCnt)
        box = np.int0(cv2.boxPoints(rect)) # vertices of the bounding rect
        center = np.int0(np.sum(box, axis=0)/4) # [centerX, centerY] dataType: int
        agent.set(center[0], center[1]) # update the position of the agent
        imageOriginal = cv2.drawContours(imageOriginal, [box], 0, (0, 255, 0), 3) # draw boundingRect on the
    return imageOriginal
```

```
class Agent():
    def __init__(self):
        self.x = 0
        self.y = 0
        self.orientation = 0
```

3. Instances of **Agent** class can be accessed via “vision.<agentName>”. Information about the agents are often used in a subthread.

# Subthread

Use a subthread when you want to apply a time-varying magnetic field with respect to the position/orientation of the agents.

```
def taskSinFieldXY(self):  
    startTime = time.time()  
    while True:  
        t = time.time() - startTime # elapsed time (sec)  
        theta = 2 * pi * self.freq * t  
        fieldX = 2 * cos(theta)  
        fieldY = 2 * sin(theta)  
        self.field.setX(fieldX)  
        self.field.setY(fieldY)  
        print('X: {}, Y: {}'.format(self.vision.agent1.x, self.vision.agent1.y))  
        if self.stopped:  
            return
```

Obtain elapsed time

Access detected objects

Stop flag

# Some useful features

- Left click on the camera image returns xy coordinate in the terminal.
- .....
- We need you to improve it!